

Meet the Perfect Liquid Handling Companion for Proteomic and Genomic Work

Increase the consistency and efficiency of genomic and proteomic liquid handling steps with the Thermo Scientific Matrix Hydra DT 96-channel automated pipettor. This compact benchtop system achieves more consistent results as compared to handheld pipetting methods.

- Transfer 96 samples at once for increased throughput
- Accommodate various protocols with a wide volume range of 0.5 to 300 ul
- Optimize pipetting performance with Thermo Scientific D.A.R.T.s* automation tips

For a limited time receive
a Hydra® DT **\$12,000**
and receive
2 Free Cases of D.A.R.T.s tips
to get you on your way to hands-free pipetting!



* D.A.R.T.s (Disposable Automation Research Tips)

Meet the Perfect Liquid Handling Companion for Proteomic and Genomic Work

For a limited time receive
a Hydra® DT **\$12,000**
and receive
2 Free Cases of D.A.R.T.s tips
to get you on your way to hands-free pipetting!



Ordering Guide and Tip Compatibility

Thermo Scientific Matrix D.A.R.T.s Tip Compatibility and Specifications

Item No.	Channels	Volume Range	Compatible Matrix D.A.R.T.s Tips	Precision	Accuracy	List	Promotion
1096-DT-100	96	0.5 µl-100 µl	0.5 µl-30 µl (5586/5587/5588)	+/-2.0% or 0.15 µl	2.0% or 0.10 µl	\$27,603	\$12,000
			XL (5506/5507/5508)	+/-2.0% or 0.15 µl	2.0% or 0.10 µl		
			1.0 µl-100 µl (5526/5527/5528)	+/-2.0% or 0.5 µl	2.0% or 0.25 µl		
1096-DT-300	96	5.0 µl-300 µl	5.0 µl-300 µl (5516/5517/5518)	+/-2.0% or 1.0 µl	2.0% or 0.75 µl	\$25,670	\$12,000
			XL (5536/5537/5538)	+/-2.0% or 1.0 µl	2.0% or 0.75 µl		
			WBXL (5546/5547/5548)	+/-2.0% or 1.0 µl	2.0% or 0.75 µl		

XL = Extended Length XLWB = Extended Length Wide Bore

Thermo Scientific Hydra DT Product Specifications

Dimensions (W x D x H)	21" / 53 cm x 14" / 35.5 cm x 10.5" / 27 cm
Volume Range	0.5 µl-100 µl; 5.0 µl-300 µl
Dispense Type	Contact
Compatible Labware	24-, 96-, 384-channel standard, low profile and deep well labware, tube racks
User Interface	On-board touch pad or ControlMate software



Contact your representative to meet your
perfect liquid handling companion and arrange
a live product demonstration in your lab today!

Offer valid in North America through April 30, 2010.

Visit www.thermo.com and enter **Hydra DT** to learn more.

© 2010 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. All other trademarks and product names are the intellectual/ industrial property of their respective owners.

North America: Tel: 800.345.0206 | email: matrix.info@thermofisher.com
Europe: Tel: +44 (0) 161 486 2110 | email: matrix.eu.info@thermofisher.com
Asia: email: matrix.ap.info@thermofisher.com

www.thermo.com/matrix

Rev B

Thermo
SCIENTIFIC