

Thermo Scientific Matrix Automation Tips are compatible with Beckman instruments using 96 format pipetting heads and are validated to ensure the highest quality product for automated pipetting.

Thermo Scientific Matrix Automation Tips



Highest Quality. All Matrix Automation Tips are molded using hydrophobic polypropylene for superior performance on 96 Beckman instruments and are validated to ensure optimal automated liquid handling.

Application Flexibility. Matrix Automation Tips are available in non-sterile, sterile and sterile filtered options, facilitating a wide range of pipetting applications on compatible instrument platforms.

Knowledgeable Support. Our team has a strong reputation for supporting both equipment and consumables, and can assist in set up and validation of Matrix Automation Tips on your 96 Beckman instrument as well as address unique supply and troubleshooting needs.

Thermo Scientific Matrix Automation Tips

Description	96 Beckman	96 Beckman
Item No.	5371, 5372, 5373	5381, 5382, 5383
Volume Range	5-250 µl	0.5-20 µl
Tip Bottom Outer Orifice	0.040" 1.016 mm	0.030" 0.762 mm
Tip Bottom Inner Orifice	0.024" 0.609 mm	0.015" 0.381 mm
Overall Length	2.029" 51.536 mm	1.154" 29.311 mm
Length Below Tip Shoulder	1.639" 41.630 mm	1.124" 28.549 mm

Ordering Guide

Item No.	Format	Description	Compatible Instruments	Quantity
5371	96	250 µl	96 Beckman	10 racks/case
5372	96	250 µl, Sterile	96 Beckman	10 racks/case
5373	96	250 µl, Sterile, Filtered	96 Beckman	10 racks/case
5381	96	20 µl	96 Beckman	10 racks/case
5382	96	20 µl, Sterile	96 Beckman	10 racks/case
5383	96	20 µl, Sterile, Filtered	96 Beckman	10 racks/case

North America: Tel: 800.345.0206 | email: matrix.info@thermofisher.com
 Europe: Tel: +44(0) 161 486 2110 | email: matrix.eu.info@thermofisher.com
 Asia: email: matrix.ap.info@thermofisher.com
www.thermo.com/matrix