

The Thermo Scientific ALTO-8 screw top decapper is designed to work exclusively with Thermo Scientific Abgene® 2D Barcoded Twist-Lock™ Storage Tubes. Increase laboratory throughput and improve ergonomics by removing and replacing screw caps on an entire column of tubes at once.

Thermo Scientific ALTO-8™ Screw Top Decapper



Improved Productivity. With The Thermo Scientific ALTO-8 screw top decapper eliminate cumbersome handwork involved with decapping and capping single tubes. Process tube racks faster by uncapping a full column of tubes at once for multichannel pipette access and compatibility with automated liquid handlers, while reducing the risk of repetitive stress injury.

Compact Design. The small, lightweight ALTO-8 decapper is easily portable from one laboratory area to another.

Superior Sample Integrity. With the 8-channel decapper design, eliminate the risk of sample contamination by preventing caps from passing over open tubes.

Knowledgeable Support. Our team has a strong reputation for supporting both equipment and consumables, and can demonstrate simple ways to improve sample management and productivity in your laboratory using the ALTO-8 decapper.

Thermo Scientific ALTO-8 Screw Top Decapper

Dimensions (W x D x H)	267 mm x 187 mm x 129 mm
Weight	9 lbs/4 kg
Power Supply	100-240 V, 47-63 Hz, 3.2 A max
Capacity	8 caps (1 column) simultaneously
Speed	15 seconds per column
Operating Humidity Range	20-85% Non-condensing
Operating Temperature Range	18 °C to 28 °C

Ordering Guide

Item No.	Description	Qty
AB-2500	ALTO-8 Screw Top Decapper	Each

Compatible Products

AB-1461	Abgene 0.5 ml Blank Twist-Lock Tubes	10 racks of 96 tubes/case
AB-1459	Abgene 0.5 ml Blank Twist-Lock Tubes w/ Caps	10 racks of 96 tubes/case
AB-1263	Abgene 1 ml Blank Twist-Lock Tubes	10 racks of 96 tubes/case
AB-1262	Abgene 1 ml Blank Twist-Lock Tubes w/ Caps	10 racks of 96 tubes/case
AB-1316	Abgene 0.5 ml Alphanumeric Twist-Lock Tubes	10 racks of 96 tubes/case
AB-1315	Abgene 0.5 ml Alphanumeric Twist-Lock Tubes w/ Caps	10 racks of 96 tubes/case
AB-1261	Abgene 1 ml Alphanumeric Twist-Lock Tubes	10 racks of 96 tubes/case
AB-1258	Abgene 1 ml Alphanumeric Twist-Lock Tubes w/ Caps	10 racks of 96 tubes/case
AB-1314	Abgene 0.5 ml 2D Barcoded Twist-Lock Tubes	10 racks of 96 tubes/case
AB-1313	Abgene 0.5 ml 2D Barcoded Twist-Lock Tubes w/ Caps	10 racks of 96 tubes/case
AB-1438	Abgene 1 ml 2D Barcoded Twist-Lock Tubes, Bulk	960 tubes/case
AB-1249	Abgene 1 ml 2D Barcoded Twist-Lock Tubes	10 racks of 96 tubes/case
AB-1248	Abgene 1 ml 2D Barcoded Twist-Lock Tubes w/ Caps	10 racks of 96 tubes/case
AB-1250	Abgene Twist-Lock Caps, Bulk	960 caps/case
AB-1317	Abgene Colored Twist-Lock Caps, Bulk, 240 each of 4 colors	960 caps/case