

Thermo Scientific Matrix CapMats and CapStrips maintain a strong, constant seal on storage tubes, microplates and deep-well blocks to ensure sample integrity. Applied with the Thermo Scientific SuperSealer, sealing is secure and consistent every time.

Thermo Scientific Matrix® CapMats and CapStrips



Seal entire plates, blocks or tube racks with Matrix CapMats or individual storage tubes with separable Matrix CapStrips.



Sample Security. Thermo Scientific Matrix CapMats, the most effective seals for long-term storage or shipping, ensure plates, blocks and tube racks remain securely sealed at temperatures down to -80°C . Simplify well identification with Alphanumeric designations in each cap well.

Sealing Flexibility. Offering the flexibility to seal and subsequently access rows of a plate, block or tube rack, or handle individual storage tubes with separable strips, Thermo Scientific Matrix CapStrips reduce the risk of cross-contamination and sample exposure while effectively maintaining a seal down to -20°C .

Highest Quality. All seals are manufactured from 100% virgin, chemically-resistant ethyl vinyl acetate.

Knowledgeable Support. Our team has a strong reputation for supporting both equipment and consumables, and can demonstrate how Matrix CapMats and CapStrips are the ideal solution for your long-term tube sealing needs.

Thermo Scientific Matrix CapMats and CapStrips

Ordering Guide

Item No.	Description	Qty
4410	CapMat for 1 ml Blocks, 96-well Plates and Snap Racks, Bulk	100/case
4411	CapMat for 1 ml Blocks, 96-well Plates and Snap Racks	10/case
4412	CapMat for 1 ml Blocks, 96-well Plates and Snap Racks, Sterile	10/case
4420	CapMat for 2 ml Blocks, Bulk	100/case
4421	CapMat for 2 ml Blocks	10/case
4422	CapMat for 2 ml Blocks, Sterile	10/case
4430	CapMat for 500 µl, 0.75 ml and 1.4 ml Storage Tubes in Latch Racks, Bulk	100/case
4431	CapMat for 500 µl, 0.75 ml and 1.4 ml Storage Tubes in Latch Racks	10/case
4432	CapMat for 500 µl, 0.75 ml and 1.4 ml Storage Tubes in Latch Racks, Sterile	10/case
4413	CapStrips, Separable for 1 ml Blocks and Storage Tubes, Bulk	240/case
4414	CapStrips, Separable for 1 ml Blocks and Storage Tubes, Bulk, Sterile	240/case
4415	CapStrips, Non-Separable for 1 ml Blocks and Storage Tubes, Bulk	240/case
4416	CapStrips, Non-Separable for 1 ml Blocks and Storage Tubes, Bulk, Sterile	240/case

*Compatibility tested with Thermo Scientific Matrix Microplates, Blocks and Storage Tubes

Compatible Products

Thermo Scientific SuperSealer