



Thermo Scientific Matrix[®] Microplate Data Sheet

Thermo Scientific Matrix 96-Well Polystyrene Flat Bottom Microplates

Technical Specifications

Item no.	Sterile	Description	Surface Properties
4915		Clear	Medium binding
4916	●	Clear	Medium binding
4922		White	Medium binding
4923		White	High binding
4924		White/Clear bottom	Medium binding
4925	●	White/Clear bottom	Medium binding
4926		Black	Medium binding
4927		Black	High binding
4928		Black/Clear bottom	Medium binding
4929	●	Black/Clear bottom	Medium binding
4930		Clear	High binding

Plate Parameters

Composition:	Polystyrene
Well bottom:	Flat
Well shape:	Cylindrical

Dimensions (mm)

Bottom length:	127.76
Top length:	123.8
Corner well center from short sides:	11.24
Bottom width:	85.48
Corner well center from long sides:	14.38
Total height:	14.4
Distance between cell centers:	9.0
Internal well dimensions – top:	6.69
Internal well dimensions – bottom:	6.58
Well depth:	10.9
Maximum volume:	370 µl
Working volume:	50 µl-300 µl

Packaging

4 sleeves of 20 plates per case; Total of 80 plates

All dimensions are in millimeters.

© 2008 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

North America: Tel: 800.345.0206 | email: matrix.info@thermofisher.com

Europe: Tel: +44 (0) 161 486 2110 | email: matrix.eu.info@thermofisher.com

Asia: email: matrix.ap.info@thermofisher.com

www.thermo.com/matrix

