

Thermo Scientific sample storage freezer racks enable compact storage of microplates, blocks and tube racks in a laboratory's cold storage unit or on the bench top.

Thermo Scientific Sample Storage Freezer Racks



Versatile Design. Thermo Scientific sample storage freezer racks are available for use with microplates, blocks and storage tube racks up to 70 mm tall. Twin positioned handles facilitate use in both chest and upright freezers and contents are kept in position during transportation with a locking rod.

Secure Storage. Maximize bench top space or enable compact storage in your refrigerator or freezer down to -80 °C with our freezer racks.

Knowledgeable Support. Our team has a strong reputation for supporting both equipment and consumables, and can determine which sample storage freezer rack is best for your laboratory's storage needs.

Trusted Resource. Researchers worldwide protect more than two billion samples inside Thermo Scientific cold storage products. With proven equipment including cryopreservation units, storage tubes, capping systems and barcode readers, you're free to concentrate on your work without worrying about your valuable samples.



Thermo Scientific Sample Storage Freezer Racks

Wall/Shelf Material	Polypropylene
Handle/Rod Material	Aluminum
Top/Base Material	Aluminum
Minimum Temperature	-80 °C

Ordering Guide

Item No.	Description	Shelves	Dimensions (W x D x H)	Compatible Products
AB-0991/G	20 mm Freezer Rack	10	105 mm x 136 mm x 240 mm	All standard height 96-, 384- and 1536-well microplates
AB-1031/G	34 mm Freezer Rack	10	105 mm x 136 mm x 383 mm	Thermo Scientific Abgene® 300 µl tube racks; Thermo Scientific Matrix® 500 µl tube racks (except 2D Barcoded ScrewTop); Thermo Scientific Matrix 384-well, deep-well blocks; and all compatible products listed above
AB-0993/G	46 mm Freezer Rack	10	105 mm x 136 mm x 505 mm	Thermo Scientific Matrix 500 µl 2D Barcoded ScrewTop tube racks; Thermo Scientific Abgene 500 µl 2D Barcoded Twist-Lock tube racks; Thermo Scientific Abgene 0.65 ml tube racks; Thermo Scientific Matrix 0.75 ml tube racks; Thermo Scientific Matrix 96-well, 1 ml and 2 ml deep-well blocks; and all compatible products listed above
AB-1122/G	58 mm Freezer Rack	10	105 mm x 136 mm x 605 mm	Thermo Scientific Abgene 1 ml 2D Barcoded Twist-Lock tube racks; Thermo Scientific Abgene 1.2 ml tube racks; Thermo Scientific Matrix 1.4 ml tube racks; and all compatible products listed above
AB-1152/G	70 mm Freezer Rack	5	105 mm x 136 mm x 383 mm	Thermo Scientific Matrix 1.0 ml 2D Barcoded ScrewTop tube racks; and all compatible products listed above

© 2008 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

North America: Tel: 800.345.0206 | email: matrix.info@thermofisher.com
 Europe: Tel: +44 (0) 161 486 2110 | email: matrix.eu.info@thermofisher.com
 Asia: email: matrix.ap.info@thermofisher.com

www.thermo.com/matrix

Thermo
SCIENTIFIC