


# Thermo Scientific Matrix® Automation Consumables




## Automation Reservoir Specifications


Item No.	Description*	Minimum Volume	Maximum Volume with Baffle	Maximum Volume without Baffle	Height of Reservoir
1064-05-8	96-Channel	5 ml	95 ml	125 ml	1.0780" / 27.40 mm
1064-15-8	96-Channel, Sterile	5 ml	95 ml	125 ml	1.0780" / 27.40 mm
1064-05-6	96-Channel, Deep Well	5 ml	160 ml	220 ml	1.6180" / 41.10 mm
1064-15-6	96-Channel, Sterile, Deep Well	5 ml	160 ml	220 ml	1.6180" / 41.10 mm
1064-05-7	384-Channel	4 ml	65 ml	95 ml	0.0696" / 17.70 mm
1064-15-7	384-Channel, Sterile	4 ml	65 ml	95 ml	0.0696" / 17.70 mm

\*packaged 20 reservoirs per case.


### Stackable



### 384 Reservoirs



### Deep Well Reservoirs



© 2008 Thermo Fisher Scientific Inc. All rights reserved.  
All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

North America: Tel: 800.345.0206 | email: [matrix.info@thermofisher.com](mailto:matrix.info@thermofisher.com)  
Europe: Tel: +44 (0) 161 486 2110 | email: [matrix.eu.info@thermofisher.com](mailto:matrix.eu.info@thermofisher.com)  
Asia: email: [matrix.ap.info@thermofisher.com](mailto:matrix.ap.info@thermofisher.com)  
[www.thermo.com/matrix](http://www.thermo.com/matrix)

P-2007-5003