

Thermo Scientific VisionMate® FlatBed 2D Barcode Reader

The Thermo Scientific VisionMate FlatBed 2D Barcode Reader decodes one to three 2D barcoded storage tube racks in 96-tube format in a single scan with ActiveX controls for database integration.



Speed and Simplicity. Operated with the VisionMate software interface provided, the Thermo Scientific VisionMate FlatBed scans up to three racks in approximately ten seconds each, making it ideal for reading a large quantity of tubes at once.

Tube Readability. The VisionMate FlatBed reader comes equipped to read all standard Thermo Scientific Matrix® 2D Barcoded Storage Tubes in 8x12 format with ease.

Seamless Integration. The VisionMate FlatBed utilizes ActiveX controls for seamless integration into your automated system or a new or existing database. The VisionMate FlatBed can also be controlled by external applications to initiate a read or acquire decoded tube data.

Data File Flexibility. With data files created as both comma-delimited linear and columnar ASCII text, choose which data format to transfer to a spreadsheet or export to your database. Fields in the VisionMate software allow for integration of linear barcodes on your tube racks into saved file names, for further sample identification.

Knowledgeable Support. Our team has a strong reputation for supporting both equipment and consumables, and can demonstrate the VisionMate FlatBed, detailing how it can streamline sample management in your laboratory.

Thermo Scientific VisionMate FlatBed 2D Barcode Reader

Speed	10-15 seconds/rack
Decoding Capacity	1-3 racks
Data Transfer Method	Data transfer to comma-delimited linear and columnar ASCII text files; ActiveX controls
Dimensions (W x D x H)	12" / 30.5 cm x 18.75" / 47.5 cm x 3" / 7.5 cm*
Power Supply	115-240 V AC, 50-60 Hz
Supported Barcode Types	Any Thermo Scientific Matrix 2D Barcoded Storage Tube in 8x12 format as well as virtually any other 2D barcode in 8x12 format
Computer Requirements	Microsoft® Windows® XP, 2 USB ports*

Ordering Guide

Item No.	Description	Qty
3111	VisionMate FlatBed 2D Barcode Reader (Includes scanner, rack template and software)	1 each

*Due to rapidly changing technology, supported dimensions and computer requirements are subject to change